

Thermo Scientific PreCon



Automated Trace Gas Pre-Concentrator

Thermo Scientific PreCon – Automat

N_2O and CH_4 are trace gases that are increasing in the atmosphere, contributing to the greenhouse effect. N_2O also affects the global ozone distribution. Monitoring the concentrations of these gases in the atmosphere has given insight into the overall global budget, but monitoring the isotope ratios of these trace gases is required for a detailed understanding of their sources.

Isotope ratio analysis of these gases has been hampered by the necessity of collecting and manually processing very large sample volumes (e.g. 70 L of air for a $^{13}C/^{12}C$ measurement of CH_4).


The Thermo Scientific PreCon allows sample sizes to be reduced by three orders of magnitude and improves sample throughput by at least one order of magnitude.

Principle of the PreCon Process

The PreCon is a fully automated trace gas pre-concentrator for the isotope characterization of N_2O and CH_4 at natural concentrations and natural isotope abundances. It attaches to any Thermo Scientific isotope ratio monitoring device like the GC Isolink or the GasBench II interfaced to a Thermo Scientific isotope ratio MS. Both CH_4 ($^{13}C/^{12}C$) and N_2O ($^{15}N/^{14}N$ and $^{18}O/^{16}O$) can be measured at natural concentrations (CH_4 : 1700 ppb; N_2O : 330 ppb) in air sample volumes <100 mL. Samples can be provided either manually using containers with valves or fully automated by using containers with septa and a Thermo Scientific GC-PAL autosampler.

For both N_2O and CH_4 the PreCon carries out sample concentration and purification (removal of H_2O , CO_2 and condensable (for CH_4) and non-condensable (for N_2O) gases), cryofocusing, and injection onto a gas chromatograph. For CH_4 , combustion in a special combustion oven with 100% combustion efficiency and a CO_2 pre-concentration step are added prior to cryofocusing. After GC separation the samples are transferred in a stream of He carrier gas into the ion source of the isotope ratio MS, and integration is carried out over the entire GC peak. Samples are bracketed by reference pulses from the automated reference gas injection system, allowing highly precise isotope ratio determination.


PreCon Schematic



$\delta^{15}N$ and $\delta^{18}O$ from N_2O

N_2O is formed as a by-product of nitrification and as an intermediate in denitrification. These processes, which operate on different N pools and have different isotopic fractionation, make distinctive contributions to the ^{15}N and ^{18}O composition of N_2O . It is thus critically important when characterizing N_2O to measure both $^{15}N/^{14}N$ and $^{18}O/^{16}O$. The PreCon analyzes N_2O as the intact molecule, so that the important $^{18}O/^{16}O$ information is preserved.

ed Trace Gas Pre-Concentrator



Trace Studies

The high precision required to carry out work at natural abundance means that the PreCon can easily handle isotopically enriched N₂O or CH₄, in which either ¹³C or ¹⁵N has been enhanced above natural levels by incorporation of tracers.

Analytical Performance


100 mL air with natural trace gas concentrations and natural isotope abundances.

EXTERNAL PRECISION (N=5)	CONC. (PPM)	δ (‰)
¹³ C/ ¹² C from CH ₄	1.7	0.5
¹⁵ N/ ¹⁴ N from N ₂ O	0.33	0.5
¹⁸ O/ ¹⁶ O from N ₂ O	0.33	1.0

PreCon Integration and Automation

The PreCon can be loaded manually or operates fully automated using the GC-PAL autosampler with a two line needle for continuous sample transfer. The GC-PAL can be equipped with a 96 x 10 mL sample tray. User-defined sample trays for higher sample volumes and specific container geometries can easily be registered and automated due to the fully object oriented structure of the GC-PAL.

All processes are controlled and can be expanded by user-definable ISL scripts in the Isodat software suite.



Installation Requirements

Instrumentation

The Thermo Scientific PreCon attaches to any Thermo Scientific *irm*-GC/MS interface, either

- GC IsoLink
- GC/GP interface
- GasBench II

The interface requires a Thermo Scientific Isotope Ratio Mass Spectrometer e. g. DELTA V series, MAT 253.

Sample Containers

Analysis of CH₄ and N₂O requires separate samples due to different operating conditions in the interface. CO is not removed prior to combustion of CH₄ so that a correction for background CO (70 ppb in air) is required for most accurate analysis.

- 100 mL flasks with gas tight manual valves
- Connectors with ¼" Cajon Ultra-Torr fittings
- 250-350 mm distance between connectors (flask length)
- Automation with the GC-PAL system requires sample containers with septa

Gases

- Ultra-pure He (≥5.0); additional purification is recommended (e.g. zeolite filled, LN₂ cooled trap)
- Calibrated CO₂ and N₂O reference gases must be provided locally.

Liquid Nitrogen

About 3 L/day

Power

230 V Single phase, 1.5 A, 154 W

Dimensions

230 mm (width) x 765 mm (depth) x 640 mm (height)
9.1 (width) x 30.2 (depth) x 25.2 inch (height)

Weight

About 23 kg (51 pound)

References

The PreCon is described and its performance is characterized in:

- ¹ W. A. Brand, PreCon: A Fully Automated Interface for the Pre-Gc Concentration of Trace Gases on Air for Isotopic Analysis, *Isotopes Environ. Health Stud.* **1995**, *31*, 277-284.
- ² A. L. Rice, et. al., High-Precision Continuous-Flow Measurements of δ¹³C and δD of Atmospheric CH₄ *Anal. Chem.* **2001**, *73*, 4104-4110.

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

Africa-Other

+27 11 570 1840 • analyze.sa@thermo.com

Australia

+61 2 8844 9500 • analyze.au@thermo.com

Austria

+43 1 333 50 34 0 • analyze.at@thermo.com

Belgium

+32 2 482 30 30 • analyze.be@thermo.com

Canada

+1 800 530 8447 • analyze.ca@thermo.com

China

+86 10 8419 3588 • analyze.cn@thermo.com

Denmark

+45 70 23 62 60 • analyze.dk@thermo.com

Europe - Other

+43 1 333 50 34 0 • analyze.emea@thermo.com

Finland / Norway / Sweden

+46 8 556 468 00 • analyze.se@thermo.com

France

+33 1 60 92 48 00 • analyze.fr@thermo.com

Germany

+49 6103 408 1014 • analyze.de@thermo.com

India

+91 22 6742 9434 • analyze.in@thermo.com

Italy

+39 02 950 591 • analyze.it@thermo.com

Japan

+81 45 453 9100 • analyze.jp@thermo.com

Latin America

+1 608 276 5659 • analyze.la@thermo.com

Middle East

+43 1 333 50 34 0 • analyze.emea@thermo.com

Netherlands

+31 76 579 55 55 • analyze.nl@thermo.com

South Africa

+27 11 570 1840 • analyze.sa@thermo.com

Spain

+34 914 845 965 • analyze.es@thermo.com

Switzerland

+41 61 716 77 00 • analyze.ch@thermo.com

UK

+44 1442 233555 • analyze.uk@thermo.com

USA

+1 800 532 4752 • analyze.us@thermo.com

www.thermo.com

Thermo Fisher Scientific (Bremen) GmbH
Management System Registered to ISO 9001:2008

©2009 Thermo Fisher Scientific Inc. All rights reserved.
All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Specifications, terms and pricing are subject to change.
Not all products are available in all countries. Please consult your local sales representative for details.

BR30188_E 10/09C

Thermo
SCIENTIFIC